virion\serion

virion\serion

YOUR  
GLOBAL  
PARTNER  
IN  
DIAGNOSTICS

SERION ELISA classic Helicobacter pylori IgA/IgG/IgM

## SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM

**Manufacturer**  
**Fabricante**  
**Κατασκευαστής**  
**Fabricante**  
**Výrobce**

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19  
D - 97076 Würzburg, Germany  
Telefon: +49 (0) 9 31 / 30 45 0  
Fax: +49 (0) 9 31 / 30 45 100

E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)

Instructions - English  
Instrucciones de empleo - Español  
Οδηγίες χρήσης - Ελληνικά  
Instruções de emprego - Português  
Pokyny - Český

(Version/Versión/Έκδοση/Versão/Verze 13.11/12-1)

# Updates

Please pay attention to the differences in comparison to the previous version.

<b><u>Current version Nr.:</u></b>	<b>V 13.11/12-1</b>
<b><u>Previous version:</u></b>	<b>V 12.10/01-1</b>
<b><u>Update in section:</u></b>	<b>General Update, 5, 7.2.1</b>

## **SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM**

### **CONTENTS**

- 1 INTENDED USE**
- 2 DIAGNOSTIC RELEVANCE**
- 3 TEST PRINCIPLE SERION ELISA *classic***
- 4 KIT COMPONENTS**
- 5 MATERIAL REQUIRED BUT NOT SUPPLIED**
- 6 STORAGE AND STABILITY**
- 7 TEST PROCEDURE SERION ELISA *classic***
  - 7.1 Evidence of Deterioration
  - 7.2 Sample Preparation and Storage
  - 7.3 Preparation of Kit Reagents
  - 7.4 Overview - Test Procedure
  - 7.5 Manual Test Procedure
  - 7.6 Automated Test Procedure
  - 7.7 Positive Control / Accuracy Control
- 8 TEST EVALUATION**
  - 8.1 Single-Point Quantification with the 4PL Method
  - 8.2 Criteria of Validity
  - 8.3 Calculation SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM
  - 8.4 Limits of Quantification
  - 8.5 Borderline Ranges
  - 8.6 Interpretation of Results
  - 8.7 Reference Range of healthy Individuals
- 9 PERFORMANCE CHARACTERISTICS**
  - 9.1 Sensitivity and Specificity
  - 9.2 Reproducibility
- 10 SAFETY MEASURES**
  - 10.1 Statements of Warning
  - 10.2 Disposal
- 11 REFERENCES**

current version No.: V 13.11/12-1  
previous version: V 12.10/01-1

EN

ES

GR

PT

CZ

# SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM

## Enzyme-immunoassay for determination of human antibodies for *in vitro* diagnostic use

SERION ELISA <i>classic</i> Helicobacter pylori IgA	Order Nr.:	ESR118A
SERION ELISA <i>classic</i> Helicobacter pylori IgG	Order Nr.:	ESR118G
SERION ELISA <i>classic</i> Helicobacter pylori IgM	Order Nr.:	ESR118M

### 1 INTENDED USE

SERION ELISA *classic* Helicobacter pylori IgA, IgG and IgM tests are quantitative and qualitative immunoassays for the detection of human antibodies in serum or plasma directed against *Helicobacter pylori*. The separate detection of individual immunoglobulin classes offers a confirmation of contact with the pathogen and the determination of the disease state.

### 2 DIAGNOSTIC RELEVANCE

*Helicobacter pylori* is a spirally formed gram negative bacterium which is highly motile due to it possessed up to seven flagellae. Microbiological identification is based on positive urea, catalase and oxidase tests as well as the absence of hippurate hydrolysis and nitrate reductase.

*Helicobacter pylori* is a very host specific organism. Other *Helicobacter* species are found in a variety of mammals such as cats, dogs, swine and mice.

Phenotypical differences between *H. pylori* isolates exclusively refer to expression/non-expression of vacuolating cytotoxin (VacA) and a second toxin which is encoded by cytotoxin-associated gene (CagA). Virulent (Type I) and non-virulent (Type II) *Helicobacter* strains can be distinguished on the basis of these phenotypical differences.

Patients with duodenal ulcers are more often infected by VacA and CagA expressing Type I *H. pylori* strains. There are, however, studies which consider that a causal relationship between *Helicobacter* infection and the associated disease syndrome of stomach cancer due to VacA and CagA producing *Helicobacter* strains is unlikely.

The person-to-person transmission mechanism for *H. pylori* is still not fully elucidated. The published literature considers the oral-oral and faecal-oral route to be likely. *H. pylori* infections are divided into two distinct groups of acute and chronic with 80 to 90 % of all gastritis cases being traceable to an *H. pylori* infection.

Diseases associated with *H. pylori* infections (ie those that can be traced back to a *H. pylori* induced gastritis) are ulcer duodeni (95 %), ulcer ventriculi (90 %) and the MALT (Mucosa Associated Lymphatic Tissue) lymphoma (60-70 %). Individuals suffering from chronic *Helicobacter pylori* infections have a six fold increased risk of developing MALT-lymphoma compared to healthy individuals.

A distinction between non-invasive and invasive methods is made in the diagnosis of *H. pylori* infections. Invasive procedures are those such as histology, urease rapid test (e.g. CLO test), microbiological techniques such as culture and molecular biologic tests (PCR). The C<sup>13</sup>-breath test and serological tests belong to the group of non-invasive methods.

European recommendations for the treatment of *Helicobacter pylori* infections were constituted by the European *Helicobacter pylori* Study Group (EHPSG) in the “Maastricht-Consensus Conference” in 1997. An eradication therapy, the “triple therapy”, was recommended for all *H. pylori* positive patients with ulcers. Furthermore all dyspeptic patients younger than 45 and without distressing symptoms should be examined with non-invasive methods such as serology and, in the case of positive results, they should also be treated.

The detection of serum antibodies can be used for therapy control after eradication treatment. In this case a period of at least six months should lapse before initiating such antibody titer control measurements to ensure that significant titer changes can be detected. Particularly in the case of IgG, the possibility of persistence for months or years should be taken into account and this puts limitations on the information gleaned from purely serological control methods. The non-invasive C<sup>13</sup> breath-test is recommended for additional control.

*Helicobacter* associated infections have to be regarded in distinct relation to the age of the patients. The prevalence per life-decade increases at 10 % in the industrialized regions of North America and Central Europe. Total prevalence of the population in this area is 40 %.

### **3 TEST PRINCIPLE SERION ELISA *classic***

The ELISA (Enzyme Linked Immunosorbent Assay) is an immunoassay, which is particularly suited to the determination of antibodies in the field of infectious serology. The reaction is based on the specific interaction of antibodies with their corresponding antigen. The test strips of the SERION ELISA *classic* microtiter plate are coated with specific antigens of the pathogen of interest. If antibodies in the patient's serum sample are present, they bind to the fixed antigen. A secondary antibody, which has been conjugated with the enzyme alkaline phosphatase, detects and binds to the immune complex. The colourless substrate p-nitrophenylphosphate is then converted into the coloured product p-nitrophenol. The signal intensity of this reaction product is proportional to the concentration of the analyte in the sample and is measured photometrically.

#### 4 KIT COMPONENTS

Test Components	Pieces / Volume
<p><b>Break apart microtiter test strips each with eight antigen coated single wells,</b> (altogether 96) <b>[MTP]</b>, 1 frame. The coating material is inactivated.</p>	12 pieces
<p><b>Standard serum (ready-to-use) [STD]</b>, Human serum in protein containing phosphate buffer; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface antigen) and anti-HCV Ab; preservative: &lt; 0.1 % sodium azide; colouring: Amaranth O.</p>	2 x 2 ml
<p><b>Negative control serum (ready-to-use) [NEG]</b>, Human serum in protein containing phosphate buffer; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface antigen) and anti-HCV Ab; preservative: &lt; 0.1 % sodium azide; colouring: Lissamin Green V.</p>	2 ml
<p><b>Anti-human IgA, IgG or IgM conjugate (ready-to-use) [APC]</b>, Anti-human IgA, IgG or IgM polyclonal antibody, conjugated to alkaline phosphatase, stabilised with protein stabilisation solution; preservative: 0.01 % methylisothiazolone, 0.01 % bromnitrodioxane.</p>	13 ml
<p><b>Washing solution concentrate (sufficient for 1000 ml) [WASH]</b>, Sodium chloride solution with Tween 20 and 30 mM Tris/HCl, pH 7,4; preservative: &lt; 0.1 % sodium azide.</p>	33,3 ml
<p><b>Dilution buffer [DILB]</b>, Protein containing phosphate buffer with Tween 20; preservative: &lt; 0.1 % sodium azide; colouring: 0.01 g/l Bromphenol blue.</p>	2 x 50 ml
<p><b>Stopping solution [STOP]</b>, 1.2 N sodium hydroxide.</p>	15 ml
<p><b>Substrate (ready-to-use) [pNPP]</b>, Para-nitrophenylphosphate in solvent free buffer; preservative: &lt; 0.1 % sodium azide (Substrate in unopened bottle may have a slightly yellow coloring, which does not reduce the quality of the product!)</p>	13 ml
<p><b>Quality control certificate with standard curve and evaluation table [INFO]</b>, (quantification of antibodies in IU/ml or U/ml).</p>	2 pages

## 5 MATERIAL REQUIRED BUT NOT SUPPLIED

- common laboratory equipment
- for the IgM detection: SERION Rf-Absorbent, order no. Z200 (20 ml)
- photometer for microtitre plates with filter, wavelength 405 nm, recommended reference wavelength 620 nm - 690 nm (e.g. 650 nm)
- incubator 37 °C
- moist chamber
- distilled water
- Click-Clips (order no. VT120)

## 6 STORAGE AND STABILITY

Reagent	Storage	Stability
Microtiter strips (coated with antigen)	unopened  after opening at 2 – 8 °C in closed aluminum bag with desiccant  <i>Strips which are not used must be stored dry in the closed aluminum bag.</i>	see expiry date;  minimum shelf-life: four weeks;  shelf-life in case of proper use and storage until expiry date
Control sera / Standard sera	after opening at 2 – 8 °C	see expiry date; 24 months as of production
Conjugate	ready-to-use solution at 2 – 8 °C <i>Avoid contamination e.g. by using sterile tips.</i>	see expiry date; 28 months as of production
Dilution buffer	Unopened  after opening at 2 – 8 °C  <i>Discard cloudy solutions.</i>	see expiry date; 36 months as of production; 24 months
Washing solution	Concentrate after opening at 2 – 8 °C working dilution at 2 – 8 °C working dilution at room temperature <i>Bottles used for the working dilution should be cleaned regularly. Discard cloudy solutions.</i>	see expiry date; 2 weeks; 1 week
Substrate	ready-to-use solution at 2 – 8 °C, stored protected from light <i>Avoid contamination e.g. by using sterile tips.</i> <i>Discard if solution turns yellow (extinction against aqua dest. &gt; 0.25 OD).</i>	see expiry date; 36 months as of production
Stopping solution	After opening at room temperature	see expiry date

## 7 TEST PROCEDURE SERION ELISA *classic*

### 7.1 Evidence of Deterioration

Only use SERION ELISA *classic* reagents when using SERION ELISA *classic* immunoassays. The components must not be exchanged for reagents of other manufacturers. Standard and control sera of SERION ELISA *classic* immunoassays are defined exclusively for the test kit to be used and must not be used in other lots. Dilution buffer, washing solution, substrate and stop solution can be used for all SERION ELISA *classic* immunoassays irrespective of the lot and the test.

There are three different conjugate concentrations for each immunoglobulin class: LOW, MEDIUM, HIGH. The classification is written on each label as follows:

e.g. IgG +            low concentrated IgG conjugate  
      IgG ++        medium concentrated IgG conjugate  
      IgG +++       high concentrated IgG conjugate

In rare cases the use of special conjugate is necessary to guarantee consistent quality of our products. Special conjugates are produced in a separate lot and do not carry the “+” sign and are not exchangeable with other conjugates.

Please pay close attention to information on labels!

Unopened, all components of the SERION ELISA *classic* tests may, if stored accordingly, be used up to the expiry dates given on the labels. Reagents may not be used after date of expiry.

Dilution or alteration of the reagents may result in a loss of sensitivity.

Avoid exposure of reagents to strong light during storage and incubation. Reagents must be tightly closed after use to avoid evaporation and contamination.

To open the aluminum bag of the microtiter plate please cut off the top of the marked side only, in order to guarantee proper reclosing. Do not use the strips if the aluminum bag is damaged or if the bag with remaining strips and desiccant was not properly reclosed.

Use aseptic techniques when removing aliquots from the reagent tubes to avoid contamination. To avoid false positive results ensure not to contact or splash the top-walls of wells while pipetting conjugate. Take care not to mix the caps of the bottles and/or vials.

Reproducibility of test results is dependent on thorough mixing of the reagents. Agitate the flasks containing control sera before use and also all samples after dilution (e.g. by using a vortex mixer).

Be sure to pipette carefully and comply with the given incubation times and temperatures. Significant time differences between pipetting the first and last well of the microtiter plate

when dispensing samples and control sera, conjugate or substrate can result in different pre-incubation times, which may influence the precision and reproducibility of the results.

Optimum results can only be achieved if the instructions are strictly followed.

The SERION ELISA *classic* immunoassay is only valid if the lot-specific validation criteria on the quality control certificate are fulfilled.

Adequate washing avoids test unspecificities. Therefore, the washing procedure should be carried out carefully. All of the flat bottom wells should be filled with equal volumes of washing buffer. At the end of the procedure ensure that the wells are free of all washing buffer in order to avoid uncontrolled dilution effects. Avoid foaming!

Take care not to damage the inscription (pathogen / antibody class) on the microtiter test strips during washing and aspiration to avoid confusion.

## 7.2 Sample Preparation and Storage

Lipaemic, hemolytic or icteric samples (serum or plasma) should only be tested with caution. Obviously contaminated samples should not be tested. Serum or plasma (EDTA, citrate, heparin) collected according to standard laboratory methods are suitable samples. Samples must not be thermally inactivated.

### 7.2.1 Dilution of Samples

Before running the test, patient samples ( $V_1$ ) must be diluted in dilution buffer ( $V_2$ ) as follows:

#### SERION ELISA *classic* Helicobacter pylori IgA/IgG

$V_1 + V_2 = 1+100$	add	10 $\mu$ l	patient's sample
	each to	1000 $\mu$ l	dilution buffer

After dilution and before pipetting into the microtiter plate the samples must be mixed thoroughly to prepare a homogenous solution.

## SERION ELISA *classic* Helicobacter pylori IgM

### Interference with rheumatoid factors

Rheumatoid factors are autoantibodies mainly of the IgM class, which preferably bind to IgG immune complexes. The presence of non-specific IgM antibodies (rheumatoid factors) can lead to false-positive results in the IgM assay. Furthermore, the possibility exists, that weak-binding pathogen-specific IgM antibodies may be displaced by stronger-binding IgG antibodies leading to a false negative IgM result. Therefore it is necessary to pretreat samples with rheumatoid factor-absorbents prior to IgM detection (SERION RF-Absorbent, Order Nr.: Z200 (20 ml/100 tests)). Rf-absorption is performed by incubation of the patient's sample in Rf-dilution buffer for 15 minutes at room temperature or over night at 4 °C. The test procedure is described in a separate instruction manual.

Before running the test, rheumatoid factor-absorbent ( $V_1$ ) must be diluted 1+4 in dilution buffer ( $V_2$ ).

$V_1 + V_2 = V_3 (1 + 4)$	add	200 $\mu$ l	Rf-absorbent
		each to 800 $\mu$ l	dilution buffer

Patient's samples ( $V_4$ ) must be diluted in this Rf-dilution buffer ( $V_3$ ):

$V_4 + V_3 = 1+100$	add	10 $\mu$ l	patient's sample
		each to 1000 $\mu$ l	Rf-dilution buffer

After dilution and before pipetting into the microtiter plate the samples must be mixed thoroughly to prepare a homogenous solution.

### 7.2.2 Sample Storage

The patient's samples should not be stored for more than 7 days at 2 – 8 °C. Extended storage is possible at  $\leq -20$  °C. Avoid repeated freezing and thawing of samples. Diluted samples can be stored at 2 – 8 °C for one week.

### 7.3 Preparation of Kit Reagents

Bring all reagents to room temperature before testing.

#### 7.3.1 Microtiter Test Strips

The microtiter test strips in frames are packed with a desiccant in an aluminum bag. Take unrequired cavities out of the frame and put them back into the aluminum bag. Close bag carefully to ensure airtight conditions.

#### 7.3.2 Control Sera / Standard Sera

Control and standard sera are ready-to-use and must not be diluted any further. For each test run - independent of the number of microtiter test strips to be used - control and standard sera must be included. The standard sera should be set up in duplicate.

Do not treat control sera with Rf-absorbent.

#### 7.3.3 Anti-human IgA, IgG or IgM AP-Conjugate (ready-to-use)

Conjugates with the same concentration and of the same immunoglobulin class are interchangeable. Avoid contamination of ready-to-use conjugates e. g. by using sterile tips.

#### 7.3.4 Washing Solution

Dilute washing buffer concentrate ( $V_1$ ) 1:30 with aqua dest. to a final volume of  $V_2$ .

Example:

Buffer concentrate ( $V_1$ )	Final volume ( $V_2$ )
33.3 ml	1000 ml
1.0 ml	30 ml

#### 7.3.5 Dilution Buffer for Samples (ready-to-use)

#### 7.3.6 Substrate (ready-to-use)

Avoid contamination of the ready-to-use substrate solution e. g. by using sterile tips.

#### 7.3.7 Stopping Solution (ready-to-use)

## 7.4 Overview - Test Procedure

### SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM quantitative

In case of IgM detection absorption of rheumatoid factor, see No. 7.2.1;  
Incubation 15 minutes at room temperature or over night at 4 °C

sample dilution<sup>1</sup>  
(patient's samples)  
1+100

Pipette diluted samples and ready-to-use control /  
standard sera into the microtest wells (100 µl)

INCUBATION 60 Min./ 37 °C  
moist chamber

WASH (4 x 300 µl [DIL] [WASH] )<sup>2</sup>

Pipette conjugate solution [APC] (100 µl)

INCUBATION 30 Min./ 37 °C  
moist chamber

WASH (4 x 300 µl [DIL] [WASH] )<sup>2</sup>

Pipette substrate solution [pNPP] (100 µl)

INCUBATION 30 Min./ 37 °C  
moist chamber

Pipette stopping solution [STOP] (100 µl)

READ EXTINCTION at 405 nm

<sup>1</sup>Special dilution buffers for the following SERION ELISA *classic* tests:  
Borrelia burgdorferi IgG, IgM, EBV EA IgG, Parvovirus B19 IgM and Hantavirus Puumala IgG, IgM

<sup>2</sup>For manual use:  
tap plate at the end of the wash procedure on paper towel.

## 7.5 Manual Test Procedure

1. Place the required number of **cavities in the frame** and prepare a protocol sheet.
2. Add each **100 µl of diluted sample or ready-to-use controls** into the appropriate wells of microtiter test strips. Spare one well for substrate blank, e.g.:

IgA/IgG/IgM quantitative	
well no.	
well A1	Substrate blank
well B1	Negative Control
well C1	Standard serum
well D1	Standard serum
well E1	Patient 1....

3. **Sample incubation** for 60 minutes (+/- 5 min) at 37 °C (+/- 1 °C) in moist chamber
4. After incubation **wash** all wells with washing solution (by automated washer or manually):
  - aspirate or shake out the incubation solution
  - fill each well with 300 µl washing solution
  - aspirate or shake out the washing buffer
  - repeat the washing procedure 3 times (altogether 4 times!)
  - dry by tapping the microtiter plate on a paper towel
5. **Addition of conjugate**  
Add 100 µl of the ready-to-use IgA/IgG/IgM conjugate to the appropriate wells (except substrate blank)
6. **Conjugate incubation** for 30 minutes (+/- 1 min)\* at 37 °C (+/- 1 °C) in moist chamber.
7. After incubation **wash** all wells with washing solution (see above)
8. **Addition of substrate**  
Add 100 µl of ready-to-use substrate solution to each well (including well for substrate blank!)
9. **Substrate incubation** for 30 minutes (+/- 1 min)\* at 37 °C (+/- 1 °C) in moist chamber.
10. **Stopping of the reaction**  
Add 100 µl stopping solution to each well, shake microtiter plate gently to mix.
11. **Read extinction**  
Read optical density (OD) within 60 minutes at 405 nm against substrate blank, reference wave length between 620 nm and 690 nm (e.g. 650 nm).

\* Please note, that under special working-conditions internal laboratory adaptations of the incubation times may be necessary.

## 7.6 Automated Test Procedure

SERION ELISA are suited for processing on automats and evaluated for use with Immunomat™ as well as with DYNEX DSX® and DS2®. The automated processing is performed analogous to manual use. Please note, that under special working-conditions internal laboratory adaptations of the incubation times may be necessary.

## 7.7 Positive Control / Accuracy Control

For the periodic verification of the test method, in order to fulfil the requirements of laboratory internal quality management systems, we recommend using SERION ELISA *controls* to determine precision and accuracy of SERION ELISA *classic* test runs. The use of SERION ELISA *controls* is described in specific instruction manuals.

## 8 TEST EVALUATION

### SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM

#### 8.1 Single-Point Quantification with the 4PL Method

Optimised assignment of extinction signals to quantitative values is guaranteed by using non-linear functions, which adjust a sigmoide curve without any further transformation to OD-values. Determination of antibody concentrations with the SERION ELISA *classic* is carried out by the 4 parameter logistic-log-model (4 PL) which is ideal for exact curve-fitting. It is based on the formula:

$$OD = A + \frac{D - A}{1 + e^{B(C - \ln \text{Conc.})}}$$

The parameters A, B, C, and D are representative for the exact shape of the curve:

- |                       |               |
|-----------------------|---------------|
| 1. lower asymptote    | ⇒ parameter A |
| 2. slope of the curve | ⇒ parameter B |
| 3. turning point      | ⇒ parameter C |
| 4. upper asymptote    | ⇒ parameter D |

For each lot the standard curve is evaluated by Institut Virion\Serion GmbH (Würzburg, Germany) in repeated test runs under optimal conditions. Time consuming and cost intensive construction of the standard curve by the user is not necessary.

For evaluation of antibody concentrations a lot specific standard curve as well as a lot specific evaluation table is included with each SERION ELISA *classic* test kit. The evaluation software SERION *evaluate* as well as the Microsoft® Excel-based software tool SERION *activity* are available on request.

To compensate for normal test variations and also for test run control a standard serum is used in each individual test run. For this control serum a reference value with a validity range is determined by the quality control of the producer. Within this range a correct quantification of antibody concentration is ensured.

## 8.2 Criteria of Validity

- The substrate blank must be < 0.25 OD.
- The negative control must produce a negative test result.
- By use of quantitative SERION ELISA *classic* tests the mean OD-value (after subtraction of the substrate blank!) of the standard serum must be within the validity range, which is given on the lot specific quality control certificate.
- The variation of OD-values of the standard serum may not be higher than 20 %.

If these criteria are not met, the test is not valid and must be repeated.

## 8.3 Calculation SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM

### 8.3.1 Non-automated Evaluation

For the SERION ELISA *classic* test evaluation a lot-specific quality control certificate with standard curve and an evaluation table is included in the test kit so that the obtained OD values may be assigned to the corresponding antibody activities. The substrate blank must be subtracted from all OD values prior to evaluation.

#### Method 1: Qualitative Evaluation

To fix the cut-off ranges multiply the mean value of the measured standard OD with the numerical data of the quality control certificate (see special case formulas), e.g.:

OD = 0.502 x MW(STD) with upper cut-off

OD = 0.352 x MW(STD) with lower cut-off

If the measured mean absorbance value of the standard serum is 0.64 OD, the range of the cut-off is in between 0.225-0.321 OD.

## Method 2:

### Continuous Determination of Antibody Activities using the Standard Curve

So called interassay variations (day to day deviations and laboratory to laboratory deviations) are compensated by multiplication of the current measured value obtained with a patient's sample with the correction factor F. This factor is calculated as follows:

$$F = \frac{\text{OD-reference value (of standard serum)}}{\text{OD-current value (of standard serum)}}$$

The procedure is necessary to adjust the current test level of the user with the lot-specific standard curve. First, daily deviations have to be corrected by calculating the correction factor F.

1. The mean of the two OD-values of the standard serum has to be calculated and checked that it is within the given validity range.
2. Calculation of the factor F: the given reference value is divided by the mean of the extinction of the standard serum:  
 $F = \text{reference value extinction STD serum} / \text{mean value extinction STD serum}.$
3. All measured values of patient's samples are multiplied by F.
4. Antibody activities in IU/ml or U/ml can be determined from the standard curve with the corrected values.

### 8.3.2 Automatic Test Evaluation with Software SERION *evaluate*

After input of the four parameters and the reference value of the standard serum, antibody activities are calculated online from processed and measured SERION ELISA *classic* test runs by the evaluation software SERION *evaluate*.

If the optical density of the standard is out of the validity range, the following message will appear.

„Standard values out of ranges in following groups: Group 1-24.” or

„Standard values differ more than 20 % in following groups: Group 1-24.”

In these cases the test run is invalid and should be repeated.

Parameters and reference value need to be changed only if there is a change of lot (evaluation table shows parameters and reference values). Correct input of the lot specific data can be checked on the basis of the standard serum activity (in IU/ml or U/ml) assigned to the standard serum. The calculated mean value of the units has to correspond to the unit value indicated on the lot specific certificate. There is an automatic correction of the measured values. In the standard version the printout displays the following:

Sample code OD-value IU/ml or U/ml Evaluation
--

## 8.4 Limits of Quantification

The limits of quantification are specified on the quality control certificate of the SERION ELISA *classic* test. The linearity of dilution within this range has been demonstrated in comprehensive evaluation studies. In case a patient sample shows a test result above the upper limit of quantification, the sample may be tested at a higher dilution. The thereby determined antibody activity must be multiplied by the additional dilution factor.

## 8.5 Borderline Ranges

The borderline ranges of the SERION ELISA *classic* Helicobacter pylori IgA/IgG/IgM tests are specified on the quality control certificates and indicate the range for borderline test results. Values obtained, when testing a patient's sample, which fall below this range indicate a negative test result; values above the borderline range are interpreted positive. In cases where the results are within the borderline range a definitive interpretation of the result is not possible. In such cases, the test should be repeated in parallel with a follow-up sample taken one to two weeks later (serum pair).

## 8.6 Interpretation of Results

Positive results in the SERION ELISA *classic* *Helicobacter pylori* IgG, IgA and IgM tests indicate the presence of specific antibodies directed against *Helicobacter pylori*. A positive serological result should always be interpreted with consideration of the patient's history and clinical picture, as a positive result does not definitely distinguish between a current disease state and antibodies due to seroprevalence rates in the general population. A negative result however indicates a seronegative status in the patient. In a case of clinical suspicion of an infection such negative results do not however rule out the possibility of infection, as the sample may have been taken too early to for antibodies to be detectable. In such cases a further sample should be taken approximately two weeks later.

*Helicobacter* serology should not be performed by focusing only on detection of IgG antibodies since the antibody constellation IgG negative/IgA positive may be found in 2 to 7 % of patients with acute *Helicobacter pylori* infections. On the other hand 15 to 35 % of patients with acute *Helicobacter pylori* infections do not produce pathogen specific IgA, making diagnosis only possible by the detection of alternative immunoglobulin classes. Therefore, both, *Helicobacter* IgG and IgA serology should be performed simultaneously. The following antibody constellations can be interpreted as follows:

**Table 1: Interpretation of various antibody constellations in *Helicobacter* serology**

IgM	IgA	IgG	interpretation
-	-	-	seronegative status, in case of clinical symptoms, control sera should be taken about 2 weeks later
+/-	-	+	seropositive status after antigen contact or active infection with lacking IgG response, seroprevalence depending on the age of patient („contagion“)
-	-	+	
+/-	+	+	suspected active infection with <i>H. pylori</i> ; typical reaction pattern in patients with (active) gastritis
+/-	+	-	suspected active infection with <i>H.pylori</i> , lacking IgG-response (2-7% of cases)

## 8.7 Reference Range of healthy Individuals

Testing of random blood donor sera, collected in the region of southern Germany, with the SERION ELISA *classic* Helicobacter pylori IgA, IgG and IgM tests resulted in the following distribution. From 47 sera tested with the SERION ELISA *classic* Helicobacter pylori IgA test, 39 (83.0 %) were negative, five sera (10.6 %) reacted positive and three sera (6.4 %) were considered borderline. Of 47 sera, 41 (87.2 %) were negative when tested with the SERION ELISA *classic* Helicobacter pylori IgG test, four (8.5 %) sera yielded a positive result and two (4.3 %) were evaluated as borderline. In addition, results from testing 47 sera in the SERION ELISA *classic* Helicobacter pylori IgM test were as follows; 35 sera (74.5 %) tested negative, six (12.8 %) positive and six (12.8 %) sera were borderline.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Sensitivity and Specificity

#### SERION ELISA *classic* Helicobacter pylori IgG:

In order to establish the performance characteristics for the SERION ELISA *classic* Helicobacter pylori IgG a total of 88 sera were examined. Of these 41 were random blood donors and 47 sera from patients with suspected *Helicobacter* infection supported by previous test results. The sera were tested in the SERION ELISA *classic* Helicobacter pylori IgG and in additional two commercially available ELISAs to provide the following performance data.

SERION ELISA <i>classic</i> Helicobacter pylori IgG	Sensitivity	Specificity
Performance data versus ELISA 1	96.6 %	94.3 %
Performance data versus ELISA 2	91.1 %	> 99 %

In addition, an external study was carried out by AFSSAPS (Agence française de sécurité sanitaire des produits de santé) in France. This consisted of the examination of a total of 92 sera. 44 of 45 samples from patients with proven Helicobacter infection were correctly identified by the SERION ELISA *classic* Helicobacter pylori test while one serum was identified as borderline. 42 of 47 serum samples from proven non-infected individuals were recorded as negative, 2 as borderline and three sera were false positive. These results support the results of the internal study. Detailed results of this external study are available at <http://www.afssaps.fr>

### SERION ELISA *classic* Helicobacter pylori IgA:

To evaluate the SERION ELISA *classic* Helicobacter pylori IgA a total of 146 serum samples were tested. Serum group 1 consisted of 54 random blood donor sera while serum group 2 contained 87 sera from gastroscopy patients whose biopsy material tested urease positive (CLO test) and were characterised as *H. pylori* positive on the basis of histological tests. In addition 5 sera from patients with antibodies against *Campylobacter jejuni* were examined for cross-reactions.

The samples were tested in the SERION ELISA *classic* Helicobacter pylori IgA in comparison to four additional commercially available ELISA tests (tests A to D). Discrepant results (i.e. any sera that did not produce the same result in at least 4 of the 5 ELISA tests) were evaluated using a commercial immunoblot test.

#### Results: Serum group 1:

Test	Correlation	Prevalence	Sensitivity	Specificity
SERION ELISA <i>classic</i>	96.0 %	20.0 %	80.0 %	100 %
ELISA A	93.0 %	16.3 %	100 %	91.7 %
ELISA B	93.9 %	20.4 %	80.0 %	97.4 %
ELISA C	84.3 %	19.6 %	20.0 %	100 %
ELISA D	95.7 %	19.1 %	88.9 %	97.4 %

A total of ten sera were determined to be positive. Eight sera were positive in the SERION ELISA *classic* Helicobacter pylori test and ELISAs B and D. ELISA A was the most sensitive test, however three sera were in the borderline region and not used in the statistical calculation. There were no significant differences in the calculations of prevalence.

#### Results: Serum group 2:

Test	Korrelation	Prevalence	Sensitivity	Specificity
SERION ELISA <i>classic</i>	91.3 %	59.4 %	85.4 %	100 %
ELISA A	82.2 %	67.1 %	100 %	45.8 %
ELISA B	88.9 %	61.1 %	81.8 %	100 %
ELISA C	86.3 %	65.0 %	78.8 %	100 %
ELISA D	89.5 %	63.2 %	83.3 %	100 %

In the SERION ELISA *classic* Helicobacter pylori IgA and ELISAs B-D, 28 sera from 87 sequentially collected sera were evaluated as definitely negative. In ELISA A 13 of these sera produced positive results. The prevalence indicates the proportion of positive results

within the serum groups. The calculated sensitivity is lower when compared to the IgG ELISA. This is partly due to the IgA immune response being primarily directed against the CagA and urease B proteins from *H.pylori* while the IgG response is directed against the complete protein complement of the organism.

The observed differences in the IgG and IgA antibody responses in addition to the fact that the urease B protein is mainly responsible for non-specific antibody binding is a probable reason for the lower correlation of the IgA determinations and, in particular, that the comparison of ELISA and immunoblot is only partly possible. This study does emphasise the importance of IgA antibody detection in the diagnosis of *H. pylori* infections and demonstrates a 2 to 5 fold higher prevalence of IgA antibodies in patients with confirmed infection (serum group 2) when compared to the normal population (serum group 1). Cross-reactions were not seen in the SERION ELISA *classic* Helicobacter pylori IgA to *Campylobacter jejuni* (five sera tested). This was also the case in ELISAs B to D whereas in ELISA A one serum was recorded as positive.

When all 146 sera are included the calculated performance characteristics for the SERION ELISA *classic* Helicobacter pylori IgA are sensitivity 84.3 % and specificity 100 % with the best correlation result of 93.5 %. Only ELISA A had a sensitivity of 100 % and this was at the cost of a low specificity of 73.8 %. ELISAs B and D are comparable to the SERION ELISA *classic* test with sensitivities of 81.5 % and 84.2 % and specificities of 98.6 % and 100 %. ELISA C had the worst sensitivity with 69.4 % (100 % specificity)

### **SERION ELISA *classic* Helicobacter pylori IgM**

To evaluate the SERION ELISA *classic* Helicobacter pylori IgM a total of 186 serum samples from adults were tested. Serum group 1 consisted of 45 random blood donor sera while serum group 2 contained 119 sera with demonstrated infection. A third serum panel, serum group 3, consisted of six sera from patients with antibodies against *Campylobacter jejuni*, seven *Chlamydia spec* seropositive, three with antibodies to *Borrelia burgdorferi*, three with *Legionella pneumophila* antibodies, one with *Brucella spec.* antibodies and ten sera with elevated rheumatoid factor activity.

The samples were tested in the SERION ELISA *classic* Helicobacter pylori IgM in comparison to two additional commercially available ELISA tests. While only two comparable ELISA tests of adequate quality and no immunoblot test were available, positive results were additionally tested in IgG and IgA ELISA tests. To determine the sensitivity and specificity, a direct comparison with the commercial test and SERION ELISA *classic* Helicobacter pylori IgM was made.

3.2 % of the serum samples (6/186) reacted positive in at least the SERION ELISA *classic* Helicobacter pylori IgM and one other ELISA test, likewise six of the sera tested were in at least one additional ELISA borderline. The international published literature indicates an IgM positive rate in the general population of 1 % to 11 %, which is supported by the results obtained in this study. In total 98 (140) sera were evaluated as negative in all ELISAs (at least in two ELISAs). This equates to 52.7 % (75.3 %). No sample reacted positive in all three ELISA tests although seven reacted borderline in at least two ELISAs (3.8 %). In total, at least 59 sera (31.7 %) reacted borderline in at least one test. Taking all 186 sera into account the SERION ELISA *classic* Helicobacter pylori IgM in comparison to

two other ELISAs delivered a sensitivity of 100 % and compared to ELISA A a specificity of 80.1 % (correlation 80.3 %) and to ELISA B a specificity of 75.5 % (correlation 76 %). The reduced sensitivities result from the fact that neither of the comparison tests recognised any sera positive, so rendering a comparison of sensitivities between these two tests impossible. If only the 45 blood donor sera are included in the calculations then the sensitivities and the correlation lies some 10 % higher. It is to note that the 45 blood donor sera were previously evaluated and designated as positive in earlier studies. In 17 of 119 patients, independent positive test results were recorded. ELISA A recorded along with the SERION ELISA *classic* Helicobacter pylori IgM a total of one positive serum and ELISA B five. IgM antibodies are frequently accompanied by IgG antibodies.

## 9.2 Reproducibility

Intraassay reproducibility was determined by testing sera of different reactivities 20 times in one test run. Interassay reproducibility was determined by testing sera of different reactivities 10 times in 10 independent assays performed on 5 different days.

$$\text{Coefficient of Variation (CV \%)} = \frac{\text{Standard deviation}}{\text{Mean value}} \times 100$$

### SERION ELISA classic Helicobacter pylori IgA:

Sample	Mean Value (OD)	Intraassay (CV %)	Mean Value (OD)	Interassay (CV %)
negative	0.499	3.4	0.467	9.9
borderline	0.627	4.6	0.580	10.0
strong positive	1.263	3.6	1.286	7.7

### SERION ELISA classic Helicobacter pylori IgG:

Sample	Mean Value (OD)	Intraassay (CV %)	Mean Value (OD)	Interassay (CV %)
borderline	0.442	4.5	0.423	13.8
strong positive	1.513	2.4	1.514	11.0
strong positive	1.615	3.1	1.599	9.8

### SERION ELISA classic Helicobacter pylori IgM:

Sample	Mean Value (OD)	Intraassay (CV %)	Mean Value (OD)	Interassay (CV %)
borderline	0.342	5.8	0.346	13.0
borderline	0.371	4.7	0.385	9.9
positive	0.658	5.5	0.661	9.2

## 10 SAFETY MEASURES

### 10.1 Statements of Warning

The SERION ELISA *classic* is designed for use by qualified personnel who are familiar with good laboratory practice.

All kit reagents and human specimens should be handled carefully, using established good laboratory practice.

- This kit contains human blood components. Although all control- and cut-off sera have been tested and found negative for anti-HIV-ab, HBs-Ag (*Hepatitis B-Virus-surface Antigen*) and anti-HCV-ab, they should be considered potentially infectious.
- Do not pipette by mouth.
- Do not smoke, eat or drink in areas in which specimens or kit reagents are handled.
- Wear disposable gloves, laboratory coat and safety glasses while handling kit reagents or specimens. Wash hands thoroughly afterwards.
- Patient's material and other potentially infectious material should be decontaminated after the test run.
- Reagents should be stored safely and be unaccessible to unauthorized access e.g. children.
- Stopping solution:  corrosive (C); causes acid burn (R34)  
Use safety glasses, gloves and laboratory coat while handling!

### 10.2 Disposal

Please observe the relevant statutory requirements!

## 11 REFERENCES

- [1] Blaser, M. J. (1998) *Helicobacter pylori* and associated diseases. *BMJ* 316, 1507-10.
- [2] Heilmann, K. L., Borchard, F. (1991) Gastritis due to spiral shaped bacteria other than *Helicobacter pylori*: clinical, histological and ultrastructural findings. *GUT* 32, 137-40.
- [3] Kosunen, T. U. *et al.* (1992) Diagnostic value of decreasing IgG, IgA and IgM antibody titres after eradication of *Helicobacter pylori*. *Lancet* 339, 893-95.
- [4] Malfertheiner, P. (1996) *Helicobacter pylori – Von der Grundlage zur Therapie* Georg Thieme Verlag Stuttgart, New York, 2. Auflage, 11-23.
- [5] Malfertheiner, P. *et al.* (1997) Current European concepts in the management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Eur. J. Gastroent. Hep.* 9, 1-2.
- [6] Malfertheiner, P. *et al.* (1998) *Helicobacter pylori* Infektionen und Ulkuserkrankungen. *Chirurg* 69, 239-48.
- [7] Malfertheiner, P. (2004) *Helicobacter pylori* Infection – ein Update 2004. *Deutsch-Medizinische Wochenschrift* 129, 1821-6.
- [8] The European *Helicobacter pylori* Study Group (EHPSG, 1997): Current European concepts in the Management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Gut* 41, 8-13.
- [9] Warren, J. R., Marshall, B. (1983) Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* 1, 1273-5.
- [10] Xiang, Z., *et al.* (1995) Analysis of expression of *cagA* and *vacA* virulence factors in 43 strains of *Helicobacter pylori* reveals that clinical isolates can be divided into two major types and that *cagA* is not necessary for expression of the vacuolating cytotoxin. *Infect Immun* 63, 94.
- [11] Yamaoka, Y. *et al.* (1999) Antibody against *Helicobacter pylori* CagA and VacA and the risk for gastric cancer. *J. Clin. Pathol.* 52, 215-8.

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- [1] Blaser, M. J. (1998) *Helicobacter pylori* and associated diseases. *BMJ* 316, 1507-10.
- [2] Heilmann, K. L., Borchard, F. (1991) Gastritis due to spiral shaped bacteria other than *Helicobacter pylori*: clinical, histological and ultrastructural findings. *GUT* 32, 137-40.
- [3] Kosunen, T. U. *et al.* (1992) Diagnostic value of decreasing IgG, IgA and IgM antibody titres after eradication of *Helicobacter pylori*. *Lancet* 339, 893-95.
- [4] Malfertheiner, P. (1996) *Helicobacter pylori – Von der Grundlage zur Therapie* Georg Thieme Verlag Stuttgart, New York, 2. Auflage, 11-23.
- [5] Malfertheiner, P. *et al.* (1997) Current European concepts in the management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Eur. J. Gastroent. Hep.* 9, 1-2.
- [6] Malfertheiner, P. *et al.* (1998) *Helicobacter pylori* Infektionen und Ulkuserkrankungen. *Chirurg* 69, 239-48.
- [7] Malfertheiner, P. (2004) *Helicobacter pylori* Infection – ein Update 2004. *Deutsch-Medicinische Wochenschrift* 129, 1821-6.
- [8] The European *Helicobacter pylori* Study Group (EHPSG, 1997): Current European concepts in the Management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Gut* 41, 8-13.
- [9] Warren, J. R., Marshall, B. (1983) Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* 1, 1273-5.
- [10] Xiang, Z., *et al.* (1995) Analysis of expression of *cagA* and *vacA* virulence factors in 43 strains of *Helicobacter pylori* reveals that clinical isolates can be divided into two major types and that *cagA* is not necessary for expression of the vacuolating cytotoxin. *Infect Immun* 63, 94.
- [11] Yamaoka, Y. *et al.* (1999) Antibody against *Helicobacter pylori* CagA and VacA and the risk for gastric cancer. *J. Clin. Pathol.* 52, 215-8.

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## 11 LITERATURA

- [1] Blaser, M. J. (1998) *Helicobacter pylori* and associated diseases. *BMJ* 316, 1507-10.
- [2] Heilmann, K. L., Borchard, F. (1991) Gastritis due to spiral shaped bacteria other than *Helicobacter pylori*: clinical, histological and ultrastructural findings. *GUT* 32, 137-40.
- [3] Kosunen, T. U. *et al.* (1992) Diagnostic value of decreasing IgG, IgA and IgM antibody titres after eradication of *Helicobacter pylori*. *Lancet* 339, 893-95.
- [4] Malfertheiner, P. (1996) *Helicobacter pylori* – Von der Grundlage zur Therapie Georg Thieme Verlag Stuttgart, New York, 2. Auflage, 11-23.
- [5] Malfertheiner, P. *et al.* (1997) Current European concepts in the management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Eur. J. Gastroent. Hep.* 9, 1-2.
- [6] Malfertheiner, P. *et al.* (1998) *Helicobacter pylori* Infektionen und Ulkuserkrankungen. *Chirurg* 69, 239-48.
- [7] Malfertheiner, P. (2004) *Helicobacter pylori* Infection – ein Update 2004. *Deutsch-Medizinische Wochenschrift* 129, 1821-6.
- [8] The European *Helicobacter pylori* Study Group (EHPSG, 1997): Current European concepts in the Management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Gut* 41, 8-13.
- [9] Warren, J. R., Marshall, B. (1983) Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* 1, 1273-5.
- [10] Xiang, Z., *et al.* (1995) Analysis of expression of *cagA* and *vacA* virulence factors in 43 strains of *Helicobacter pylori* reveals that clinical isolates can be divided into two major types and that *cagA* is not necessary for expression of the vacuolating cytotoxin. *Infect Immun* 63, 94.
- [11] Yamaoka, Y. *et al.* (1999) Antibody against *Helicobacter pylori* CagA and VacA and the risk for gastric cancer. *J. Clin. Pathol.* 52, 215-8.

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## 11 ODKAZY

- [1] Blaser, M. J. (1998) *Helicobacter pylori* and associated diseases. *BMJ* 316, 1507-10.
- [2] Heilmann, K. L., Borchard, F. (1991) Gastritis due to spiral shaped bacteria other than *Helicobacter pylori*: clinical, histological and ultrasstructural findings. *GUT* 32, 137-40.
- [3] Kosunen, T. U. *et al.* (1992) Diagnostic value of decreasing IgG, IgA and IgM antibody titres after eradication of *Helicobacter pylori*. *Lancet* 339, 893-95.
- [4] Malfertheiner, P. (1996) *Helicobacter pylori* – Von der Grundlage zur Therapie Georg Thieme Verlag Stuttgart, New York, 2. Auflage, 11-23.
- [5] Malfertheiner, P. *et al.* (1997) Current European concepts in the management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Eur. J. Gastroent. Hep.* 9, 1-2.
- [6] Malfertheiner, P. *et al.* (1998) *Helicobacter pylori* Infektionen und Ulkuserkrankungen. *Chirurg* 69, 239-48.
- [7] Malfertheiner, P. (2004) *Helicobacter pylori* Infection – ein Update 2004. *Deutsch-Medizinische Wochenschrift* 129, 1821-6.
- [8] The European *Helicobacter pylori* Study Group (EHPSG, 1997): Current European concepts in the Management of *Helicobacter pylori* infection. The Maastricht Consensus Report. *Gut* 41, 8-13.
- [9] Warren, J. R., Marshall, B. (1983) Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* 1, 1273-5.
- [10] Xiang, Z., *et al.* (1995) Analysis of expression of *cagA* and *vacA* virulence factors in 43 strains of *Helicobacter pylori* reveals that clinical isolates can be divided into two major types and that *cagA* is not necessary for expression of the vacuolating cytotoxin. *Infect Immun* 63, 94.
- [11] Yamaoka, Y. *et al.* (1999) Antibody against *Helicobacter pylori* CagA and VacA and the risk for gastric cancer. *J. Clin. Pathol.* 52, 215-8.

**SERION ELISA classic**

(V 11/12-1)

Symbole auf den Etiketten/ symbols on labels/ symboles et étiquettes/ simboli sulle etichette/ символы на этикетках/símbolos sobre las etiquetas/ σύμβολα στις ετικέτες/ símbolos nos rótulos / Symboly na štítcích / symboler på etiketter/ symboler på etiketterna/ Symbole na etykietach/ symboly na označení/ Simboli na oznakah/ symbol på etiketter

Hersteller/ Manufacturer/ Fabricant/ Produttore/Производитель/ Fabricante/ Κατασκευαστής/ Fabricante/ Výrobce/ Fremstiller/ Tillverkare/ Producent/ Výrobca/ Izdelovalec/ Produsent

Ausreichend für 96 Tests/ sufficient for 96 tests/ suffisant pour 96 tests/ sufficiente per 96 test/ достаточно для 96 тестов / suficiente para 96 pruebas/ επαρκεί για 96 δοκιμασίες/ suficiente para 96 ensaios/ stačí na 96 testů/ nok til 96 test/ tillräckligt för 96 tester/ Wystarcza na 96 testów/ postačuje na 96 testov/ Zadostuje za 96 testov/ Tilstrekkelig til 96 tester

Charge/ lot/ lot / lotto/ lote/ παρτίδα/ lote/ šarže/ lot/ lot/ seria/ šarža/ serija/ lot /lot


Referenz oder Bestellnummer/ reference or order number/ numéro de référence ou de commande/ numero di riferimento o ordinazione/ ссылка или номер для заказа / referencia o número de pedido/ Αριθμός αναφοράς ή παραγγελίας/ referència ou número para encomenda/ reference nebo číslo objednávky/ reference eller bestillingsnummer/ referens eller beställningsnummer/ Numer referencyjny lub numer zamówienia/ referenčné číslo alebo číslo objednávky/ referenčna ali kataloška številka/ Referanse eller ordrenummer

Lagern zwischen 2 und 8 Grad Celsius/ store between 2 and 8 degree celsius/ entre 2 et 8 degré celsius/ conservare a temperatura compresa tra 2 e 8 gradi centigradi/ хранить при температуре от 2 до 8 градусов цельсия / conservar entre 2 y 8 grados celsius/ Φύλαξη μεταξύ 2 και 8 βαθμούς Κελσίου/ Armazenar entre 2º e 8º Celsius/ uchovávejte při teplotě 2 až 8 °C/ opbevares mellem 2 og 8 grader celsius/ förvara vid 2 till 8 grader Celsius/ Przechowywać w temp. pomiędzy 2 a 8 stopni Celsjusza/ skladovať pri teplote 2 až 8 stupňov Celzia/ Shranjujte pri temperaturi od 2 do 8 C/ Oppbevares mellom 2 og 8 grader Celsius

CE-Markierung bei Erfüllung der IVD Richtlinie 98/79 EG/ CE marking according to IVD guideline 98/79 EC/ Étiquetage CE selon les directives DIV/ marcatura CE in conformità alla direttiva IVD 98/79 EC/ маркировка CE согласно директивам IVD 98/79 /marca CE según la directiva IVD 98/79 CE/ Σήμανση CE σύμφωνα με την οδηγία IVD 98/79 EE/ Marcação CE de acordo com a Directiva 98/79/ značení CE podle směrnice IVD 98/79/ES/ CE-mærkning iht. IVD-retningslinje 98/79/EF/ CE-mærkning enligt riktlinjerna för IVD i direktiv 98/79/EC/ Oznakowanie CE zgodne z wytycznymi dot. diagnostyki in vitro 98/79 EC/ označenie CE podľa smernice IVD 98/79/ES/ oznaka CE, skladna s smernico IVD 98/79/ES/ CE-merking i henhold til IVD-retningslinjer 98/79/EØF

CE-Markierung bei Erfüllung der IVD Richtlinie 98/79 EG gemäß Anhang II, Liste B/ CE marking according to IVD guideline 98/79 EC according to annex II, list B/ Étiquetage CE selon les directives DIV 98/79 CE selon l'annexe II, liste B/ marcatura CE in conformità alla direttiva IVD 98/79 EC secondo l'allegato II, elenco B/ маркировка CE согласно директивам IVD 98/79, приложение II, список B / marca CE según la directiva IVD 98/79 CE de acuerdo con el anexo II, lista B/ Σήμανση CE σύμφωνα με την οδηγία IVD 98/79 EE, σύμφωνα με το παράρτημα II, κατάλογο B/ Marcação CE de acordo com a Directiva 98/79/ CE relativo aos dispositivos médicos de diagnóstico *in vitro*, segundo a lista B do anexo II/ značení CE podle směrnice IVD 98/79/ ES podle příloh II, seznamu B/ CE-mærkning iht. IVD-retningslinje 98/79 /EF iflg. annekts II, liste B/ CE-mærkning enligt riktlinjerna för IVD i direktiv 98/79/EC, bilaga II, lista B/ Oznakowanie CE zgodne z wytycznymi dot. diagnostyki in vitro 98/79 EC, zgodnie z aneksem II, lista B/ označenie CE podľa smernice IVD 98/79 ES v znení dodatku II, zoznam B/ oznaka CE, skladna s smernico IVD 98/79/ES in seznamom B v Dodatku II/ CE-merking i henhold til IVD-retningslinjer 98/79/EØF, tillegg II, liste B

Verfallsdatum/ expiry date/ date d'expiration/ data di scadenza/ срок годности до /fecha de caducidad/ ημερομηνία λήξης/ data de validade/ datum exspirace/ udløbsdato/ förfallodatum/ data upływu ważności/ dátum expirácie/ datum izteka roka uporabnosti/ utløpsdatp

Mikrotiterplatte (brechbare Streifen)/ microtiter plate (breakable strips)/ plaque de microtitration (bandelettes détachables)/ piastra per microtitolazione (strisce separabili)/ микротитровальная панель (отрывные стрипы) /placa de microtitulación (tiras rompibles)/ Πλάκα μικροτιτλοποίησης (αποσπώμενες ταινίες)/ placa de microtitulação (tiras quebráveis)/ mikrotitrační deska (rozlomitelné proužky)/ mikrotiterplade (afbrækkelige strimler)/ mikrotiterplatta (brytbara strips)/ Płytką mikrotitracyjną (paski do odrywania)/ mikrotitračná platnička (rozlomitelné prúžky)/ vsebnik za mikrotitriranje (z razdelki, ki jih je mogoče odlomiti)/ Mikrotiterplatte (avbrytbare strips)

Antigen/ antigen/ Antigène/ antigene/ антиген /antígeno/ αντιγόνο/ antigénio/ antigen/ antigen/ antigen/ Antygen/ antigén/ antigen/ Antigen

- AK** Antikörper/ antibodies/ Anticorps/ anticorpi/ антитела / anticuerpos/ αντίσωμα/ anticorpos/ protilátky/ antistoffer/ antikroppar/ Przeciwciała/ protilátky/ protitelesa/ Antistoffer
- CAG** Kontrollantigen/ control antigen/ antigène de contrôle/ antigene di controllo/ контрольный антиген /antígeno de control/ αντιγόνο ελέγχου/ antígeno de controle/ kontrolní antigen/ kontrolantigen/ kontrollantigen/ antygen kontrolny/ kontrolný antigén/ kontrolni antigen/ kontrollantigen
- STD** Standardserum/ standard serum/ Sérum standard/ siero standard/ стандартная сыворотка /suero patrón/ πρότυπος ορός/ soro padrão/ standardní sérum/ standardserum/ standardserum/ Surowica standardowa/ štandardné sérum/ standardni serum/ Standardserum
- POS** Positivkontrolle/ positive control/ Contrôle positif/ controllo positivo/ положительные контроли /control positivo/ θετικός έλεγχος/ controlo positivo/ pozitivní kontrola/ positiv kontrol/ positiv kontroll/ Kontrola pozytywna/ pozitivna kontrola/ pozitivna kontrola/ Positiv kontroll
- C/O** Grenzwertiges Serum/ cut-off serum/ Sérum seuil/ siero cut-off/ сомнительные сыворотки (пограничные)/suero de corte/ οριακός ορός (cut-off)/ soro *cut-off*/ cut-off sérum/ cutoff-serum/ cutoff-serum/ Surowica „cut-off“/ sérum na určenie hraničnej hodnoty/ mejni serum/ Stopps serum
- NEG** Negativkontrolle/ negative control/ Contrôle négatif/ controllo negativo/ отрицательные контроли /control negativo/ αρνητικός έλεγχος/ controlo negativo/ negativní kontrola/ negativ kontrol/ negativ kontroll/ Kontrola negatywna/ negatívna kontrola/ negativna kontrola/ Negativ kontroll
- APC** Alkalisches Phosphatase Konjugat antihuman/ alkaline phosphatase conjugate anti-human/ conjugué phosphatase alcaline anti-humain/ coniugato con fosfatasi alcalina anti-umano/ античеловеческий щелочной конъюгат фосфатазы / conjugado anti humano de fosfatasa alcalina/ Σύζευξη αλκαλικής φωσφατάσης/ conjugado anti-humano com fosfatase alcalina/ konjugát alkalické fosfatázy anti-humánní/ alkalisk phosphatase konjugat antihumant/ antihumant alkaliskt fosfatas-konjugat/ Anty-ludzki koniugat fosfatazy alkalicznej/ konjugát antihumánnejalkalickej fosfatázy/ konjugat alkalne fosfataze, antihumani/ Alkalisk fosfatase-konjugat, anti-humant
- +** niedrig-konzentriertes Konjugat/ conjugate with low concentration/ conjugué à faible concentration/ coniugato a concentrazione bassa/ конъюгат низкой концентрации /conjugado con concentraci6n baja/ Σύζευξη χαμηλής συγκέντρωσης/ conjugado de baixa concentraç6o/ konjugát s nízkou koncentrací/ konjugat med lav koncentration/ konjugat med låg koncentration/ koniugat o niskim stężeniu/ konjugát so strednou koncentráciou/ konjugat z majhno koncentracijo/ Konjugat med lav konsentrasjon
- ++** mittel-konzentriertes Konjugat/ conjugate with medium concentration/ conjugué à concentration moyenne/ coniugato a concentrazione media/ конъюгат средней концентрации /conjugado con concentraci6n media/ Σύζευξη μέτριας συγκέντρωσης/ conjugado de concentraç6o intermédia/ konjugát se střední koncentrací/ konjugat med medium koncentration/ konjugat med medelhög koncentration/ koniugat o średnim stężeniu/ konjugát so strednou koncentráciou/ konjugat s srednj6 koncent/ Konjugat med middels konsentrasjon
- +++** hoch-konzentriertes Konjugat/ conjugate with high concentration/ conjugué à concentration élevée/ coniugato a concentrazione alta/ высококонцентрированный конъюгат/ conjugado con concentraci6n alta/ σύζευξη υψηλής συγκέντρωσης/ conjugado de elevada concentraç6o/ konjugát s vysokou koncentrací/ konjugat med høj koncentration/ konjugat med hög koncentration/ koniugat o wysokim stężeniu/ konjugát s vysokou koncentráciou/ konjugat z veliko koncentracijo/ Konjugat med høy konsentrasjon
- RF** Rheumafaktor-Absorbens (Rf-Absorbens)/ rheumatoid factor absorbent (rf-absorbent)/ absorbant de facteur rhumatoïde (rf-absorbant)/ adsorbente del fattore reumatoide (adsorbente Rf)/ абсорбент ревматоидного фактора (Rf-абсорбент) /adsorbente de factor reumatoide (material absorbente de Rf)/ Απορροφητής ρευματοειδούς παράγοντα (απορροφητής Rf)/ absorbente de factor reumatóide (absorbente de Fr)/ absorbent revmatoidního faktoru (rf-absorbent)/ reumafaktor- absorptionsmiddel (rf-absorptionsmiddel)/ reumafaktor-absorptionsmedel (rf-absorptionsmedel)/ Absorbent czynnika reumatoidalnego (absorbent RF)/ absorbent reumatoidného faktora (absorbent rf)/ absorbent revmatoidnega faktorja (absorbent RF)/ Revmatoid faktor-absorbent (rf-absorbent)
- DILB** Verdünnungspuffer für Serum/ dilution buffer for sera/ sérum pour le tampon de dilution/ tampone di diluizione per sieri / разбавляющий буфер для сыворотки / solución amortiguadora para los sueros/ ρυθμιστικό διάλυμα αραιώσης για ορούς/ tampão de diluição para soro/ ředici pufr pro séra/ fortyndingsbuffer til sera/ spädningsbuffert för serum/ bufor rozcieńczający do surowic / pufo na riedenie sér/ pufo za redčenje seruma/ Fortynningsbuffer til serum
- DILBS1**
- DILBS2**
- DILBS3**

- WASH** Waschlösungskonzentrat/ washing solution concentrate/ concentré de solution de lavage / soluzione di lavaggio concentrata / промывочный концентрат /concentrado de solución de lavado/ συμπύκνωμα έκπλυσης/ concentrado de solução de lavagem/ koncentrát promývacího roztoku/ vaskeopløsningskoncentrat/ tvättlösningkoncentrat/ Stężony roztwór do płukania/ koncentrát premývacieho roztoku/ koncentrat za raztopino za izpiranje/ Vaskeløsningskonsentrat
- pNPP** pNPP Substrat/ pNPP substrate/ substrat Pnpp/ substrato pNPP/ pNPP субстрат / sustrato pNPP/ Υπόστρωμα pNPP/ substrato pNPP/ pNPP substrát/ pNPP-substrat/ pNPP-substrat/ Substrat pNPP/ substrát pNPP/ substrat pNPP/ pNPP-substrat
- STOP** Stoppløsning/ stopping solution/ solution d'arrêt/ soluzione di arresto/стоп-раствор/ solución de parada/ διάλυμα διακοπής/ solução de paragem/ zastavovací roztok/ stopopløsning/ stoppløsning/ roztwór zatrzymujący reakcję/ ukončovací roztok/ raztopina za ustavitve reakcije/ stoppeløsning
-  ätzend/ corrosive/ corrosif/ corrosivo/ едкий /corrosivo/ διαβρωτικό/ corrosivo/ žíravý/ ætsende/ frätande/ czynnik korozyjny/ korozivne/ jedko/ etsende
-  ätzend/ corrosive/ corrosif/ corrosivo/ едкий /corrosivo/ διαβρωτικό/ corrosivo/ žíravý/ ætsende/ frätande/ czynnik korozyjny/ korozivne/ jedko/ etsende
- INFO** Gebrauchsanweisung, Zertifikat (Standardkurve und Auswertetabelle), CD/ instructions, certificate (standard curve and evaluation table), CD/ instructions, certificat (courbe de référence et tableau d'évaluation), CD/ istruzioni per l'uso, certificato (curva standard e tabella interpretativa), CD/ Инструкция по применению, сертификат (стандартная кривая и таблица для оценки), компактний диск /instrucciones, certificado (curva patrón y tabla de evaluación), CD/ Οδηγίες χρήσης, Πιστοποιητικό (πρότυπη καμπύλη και πίνακας υπολογισμού), CD/ instruções, certificado (curva padrão e tabela de avaliação), CD/ (standardní křivka a vyhodnocovací tabulka), CD/ brugsanvisning, certifikat (standardkurve og evalueringstabel), CD/ instruktioner, certifikat (standardkurva och utvärderingstabel), CD/ Instrukcje, certyfikat (krzywa standardowa i tabela do określania wyników/ CD/ pokyny, certifikát (šstandardná krivka a hodnotiacia tabuľka), disk CD/ navodila, certifikat (standardna krivulja in ocenjevalna tabela), CD/ Instruksjoner, certifikat (standardkurve og evalueringstabel), CD
- RTU** gebrauchsfertig/ ready-to-use/ prêt à l'emploi/ pronto per l'uso/ готовый к использованию /listo para usar/ έτοιμο προς χρήση/ pronto a utilizar/ připravený k použití/ klar til brug/ bruksfærdig/ gotowy do użycia/ pripravené na použitie/ pripravljen za uporabo/ klar til bruk
- CONC** Konzentrat/ concentrate/ concentré/ concentrato/ концентрат / concentrado/ Συμπύκνωμα/ concentrado/ koncentrát/ koncentrat/ koncentrat/ Konzentrat/ koncentrát/ koncentrat/ Konsentrat
- DIL** verdünnen oder lösen in/ dilute or dissolve in/ diluez ou dissoudre dans/ diluire o sciogliere in/ разбавить или растворить в /diluir o disolver en/ αραιώση ή διάλυση σε/ diluir ou dissolver em/ naředte nebo rozpustte v/ fortynd eller opløs i/ späd eller lös i/ Rozcieńczyć lub rozpuścić w/ rozriediť alebo rozpustiť v/ razredčite ali raztopite v/ Fortynnes eller løses opp i
- AQUA** destilliertes Wasser/ aqua detillata/ eau distillée/ acqua distillata/ дистиллированная вода /agua destilada/ αποσταγμένο νερό/ água destilada/ destilovaná voda/ destilleret vand/ destillerat vatten/ woda destylowana/ destilovaná voda/ destilirana voda/ Destillert vann
- IVD** In-vitro Diagnostik Anwendung/ in-vitro diagnostic use/ utilisation en diagnostic in-vitro/ uso diagnostico in vitro/ использование в диагностике ин-витро /uso diagnóstico in-vitro/ Διάγνωση, χρήση in-vitro/ para diagnóstico in vitro/ diagnostické použitie in-vitro/ til in-vitro diagnostik/ in vitro-diagnostisk användning/ do diagnostyki in vitro/ diagnostické použitie in-vitro/ uporaba pri diagnostiki in vitro/ In vitro-diagnostisk bruk

**SERION ELISA classic**

102	Masern Virus / Measles Virus / Rougeole	1262	Parainfluenza Virus 2
103	Mumps Virus / Parotitis virus / Oreillons	1263	Parainfluenza Virus 3
104	Varicella-Zoster Virus (VZV)	127	Mycoplasma pneumoniae
105	Herpes simplex Virus 1/2	128	Adenovirus
1051	Herpes simplex Virus 1	129	Röteln Virus / Rubella virus / virus de rubéole
1052	Herpes simplex Virus 2	130	Diphtherie / Diphtheria
106	Legionella pneumophila 1-7	1311	Coxiella burnetii (Q-Fieber) Phase 1 / Coxiella burnetii (Q-fever) phase 1
107	Echinococcus	1312	Coxiella burnetii (Q-Fieber) Phase 2 / Coxiella burnetii (Q-fever) phase 2
108	Tetanus	132	Aspergillus fumigatus
109	Cytomegalovirus	133	Enterovirus
110	Toxoplasma gondii	134	Coxsackievirus
112	FSME Virus / TBE Virus	135	Echovirus
113	Resp. Syncytial Virus (RSV)	1361	Epstein-Barr Virus VCA
114	Dengue Virus	1362	Epstein-Barr Virus EBNA 1
116	Brucella	1363	Epstein-Barr Virus Early Antigen
117	Candida albicans	137	Chlamydia
118	Helicobacter pylori	1371	Chlamydia pneumoniae
120	Bordetella pertussis	1372	Chlamydia trachomatis
1201	Bordetella pertussis Toxin	138	Yersinia
121	Borrelia burgdorferi	139	Campylobacter jejuni
122	Parvovirus B19	140	Bacillus anthracis
1231	Influenza A Virus	142	Francisella tularensis
1232	Influenza B Virus	143	Tyrosinase
125	Leptospira	144	Thyroglobulin
126	Parainfluenza Virus 1, 2, 3	145	Hantavirus Puumala
1261	Parainfluenza Virus 1		