

CRITERIA FOR SELECTING A MATTRESS



	FOAM MATTRESS	AIR MATTRESS WITHOUT MOTOR	AIR MATTRESS WITH MOTOR
	Optimal 5zon	CuroCell AREA Zone	CuroCell A4 CX15
Patient is not at risk for pressure ulcer	Recommended		
Patient is at low risk for pressure ulcer	Recommended	Recommended	
Patient is at high risk for pressure ulcer		Recommended	Recommended
Patient has pressure ulcer stage 1 or 2		Recommended	Recommended
Patient has pressure ulcer stage 3 or 4 (also suspected deep tissue injury and non stagable)			Recommended
Patient has unstable spinal injury	Recommended	Recommended	Do not use
Patient's weight is up to 180 kg	Recommended	Recommended	Recommended
Patient's weight is up to 220 kg		Recommended	Recommended
Patient's height is under 120 cm	Recommended	Recommended	

Remember that the patient needs to be moved and positioned regularly, no matter what kind of mattress he is using

Remember to unload the pressure from the heels by raising the heels from the mattress